


GST-Tag Mouse Monoclonal Antibody

Cat. No: DTG01

DESCRIPTION

Product Name	GST-Tag Mouse Monoclonal Antibody
Host Species	Mouse
Purification	Protein A affinity column
Applications	WB/ELISA
Specificity	Anti-GST Monoclonal Antibody (Mouse) recognizes GST tagged proteins.
Immunogen	GST tag
Isotope	IgG1
Conjugate	Unconjugated
Concentration	1.0mg/ml, monoclonal
Storage	For up to 12 months at -80°C. After received, the antibody should be divided into the right amount. Avoid repeated freezing and thawing cycles.
Storage buffer	PBS, pH7.4
Dilution	WB: 1 : 500-1 : 1000 ELISA: 0.1-1.0 µg/mL

APPLICATION



Lane 1: GST Tagged Protein (0.50 µg)
 Lane 2: GST Tagged Protein (0.10 µg)
 M: Western Blot Marker
 Using Anti- GST antibody