

Rabbit Anti Streptavidin Antibody, pAb

Cat.No: DTA0231

DESCRIPTION

Name	Rabbit Anti Streptavidin Antibody, pAb
Immunogen	Recombinant streptavidin
Host Species	Rabbit
Specificity	This antibody is specific to streptavidin.
Form	Liquid
Concentration	1 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 50% glycerol and 0.03% Proclin 300.
Purification	Antigen affinity purification
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	0.05-0.5μg/ml
Western Blot	0.1-1μg/ml
Other applications	User-optimized

Examples




Figure 1. Detection of streptavidin fusion protein were analyzed by Western blot using Rabbit Anti Streptavidin Antibody, pAb (Cat.DTA0231 , at 1 $\mu\text{g/mL}$ concentration).

The signal was developed with Goat Anti-Rabbit IgG(H&L)[HRP], pAb(Cat.DTA0252).

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!