


Interleukin-21 Protein

Cat.No:DTA0686

DESCRIPTION

Name	Interleukin-21 protein
Purification	The soluble target protein Interleukin-21 was derived from secreted expression and was purified by affinity chromatography
Purity	>85% by SDS-PAGE
Molecular weight	15.7 kDa
Expression System	HEK293
Species	Dog
Physical Appearance	Liquid
Buffer	1×PBS, 5% Glycerol, pH7.4
Concentration	0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)

SDS-PAGE & Western Blot



Lane 1: BSA (1.00 μ g)
Lane 2: IL21 (犬) protein 0.20 μ g
M₁: SDS-PAGE Marker
M₂: Western Blot Marker
Using Anti-His antibody

Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only