



NACHT, LRR and PYD domains-containing protein 3

Cat.No:DTA0699

DESCRIPTION

| | |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Name | NACHT, LRR and PYD domains-containing protein 3 |
| Purification | The soluble target NACHT, LRR and PYD domains-containing protein 3 was derived from secreted expression and was purified by affinity chromatography |
| Purity | >85% by SDS-PAGE |
| Molecular weight | 120.3kDa |
| Expression System | HEK293 |
| Species | Human |
| Physical Appearance | Liquid |
| Buffer | 1×PBS, 5% Glycerol, pH7.4 |
| Concentration | 0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford |
| Stability & Storage | For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws) |

SDS-PAGE & Western Blot



Lane 1: BSA (1.00 μg)
 Lane 2: NLRP3 protein (1.00μg)
 M₁: SDS-PAGE Marker
 M₂: Western Blot Marker
 Using Anti-His antibody

Thank you for choosing Detai Biological. We are dedicated to serving you!
 For Research Use Only