


TNF receptor superfamily member 1A Protein

Cat.No:DTA0700

DESCRIPTION

Name	TNF receptor superfamily member 1A protein
Purification	The soluble target protein TNF receptor superfamily member 1A was derived from secreted expression and was purified by affinity chromatography
Purity	>85% by SDS-PAGE
Molecular weight	19.1 kDa
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, 5% Glycerol, pH7.4
Concentration	0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)

Restriction Digestion Map



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only