


# Insulin Protein

Cat.No:DTA0701

## DESCRIPTION

<b>Name</b>	<b>Insulin protein</b>
<b>Purification</b>	The soluble target protein Insulin was derived from secreted expression and was purified by affinity chromatography
<b>Purity</b>	>85% by SDS-PAGE
<b>Molecular weight</b>	0.95 kDa
<b>Expression System</b>	HEK293
<b>Species</b>	Human
<b>Physical Appearance</b>	Liquid
<b>Buffer</b>	1×PBS, 5% Glycerol, pH7.4
<b>Concentration</b>	0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
<b>Stability &amp; Storage</b>	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)

## SDS-PAGE



**Thank you for choosing Detai Biological. We are dedicated to serving you!**  
**For Research Use Only**