


Cytochrome C Protein

Cat.No:DTA0702

DESCRIPTION

Name	Cytochrome C protein
Purification	The soluble target protein Cytochrome C was derived from secreted expression and was purified by affinity chromatography
Purity	>85% by SDS-PAGE
Molecular weight	11.7 kDa
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, 5% Glycerol, pH7.4
Concentration	0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)

SDS-PAGE



Lane 1: Cytochrome C Protein

M1: SDS-PAGE Protein Marker

Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only