


Coagulation factor II, thrombin Protein

Cat.No:DTA0705

DESCRIPTION

Name	Coagulation factor II, thrombin protein
Purification	The soluble target protein Coagulation factor II, thrombin was derived from secreted expression and was purified by affinity chromatography
Purity	>85% by SDS-PAGE
Molecular weight	20.2 kDa
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, 5% Glycerol, pH7.4
Concentration	0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)

SDS-PAGE & Western Blot



Lane 1: Custom Protein
 M₁: SDS-PAGE Protein Marker
 M₂: Western-Blot Protein Marker
 (Using Anti-His antibody)

Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only