

Mouse Anti SUMO Tag Antibody, mAb

Cat.No: DTA0714


Name	Mouse Anti SUMO Tag Antibody, mAb
Immunogen	Recombinant SUMO Tag protein
Host Species	Mouse
Specificity	This antibody recognizes SUMO tagged fusion proteins.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 0.1% BSA, 50% glycerol and 0.03% Proclin 300.
Purification	Purified by Protein A
Clonality	Monoclonal
Isotype	Mouse IgG2a, _κ
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

DESCRIPTION

Applications

Application	Recommended Usage
ELISA	0.1-1 μg/ml
Western Blot	0.1-1μg/ml
Other applications	User-optimized

Examples



Western blot analysis of SUMO-tagged fusion protein using Mouse Anti SUMO Tag Antibody, mAb (1 μg/ml). The signal was developed with Goat Anti-Mouse IgG(H&L)[HRP], pAb.

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!