


Granulocyte-macrophage colony-stimulating factor

Protein

Cat.No:DTD0037

DESCRIPTION

| | |
|--------------------------------|---|
| Name | Granulocyte-macrophage colony-stimulating factor protein |
| Purification | The soluble target protein Granulocyte-macrophage colony-stimulating factor was derived from secreted expression and obtained by affinity chromatography purification |
| Purity | >90% by SDS-PAGE |
| Molecular weight | 15 kDa |
| Expression System | HEK293 |
| Species | Pig |
| Physical Appearance | Liquid |
| Buffer | 1×PBS, pH7.4 |
| Concentration | At 1.00 mg/ml, BSA was used as the standard and the protein concentration was determined by Bradford assay |
| Stability & Storage | For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws) |



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only