


# Desmoglein-3(human) Protein

Cat.No:DTD0042

## DESCRIPTION

|                                |   |
|--------------------------------|---|
| <b>Name</b>                    | <b>Desmoglein-3(human) protein</b>  |
| <b>Purification</b>            | The soluble target protein Desmoglein-3(human) was derived from secreted expression and was purified by affinity chromatography |
| <b>Purity</b>                  | >85% by SDS-PAGE  |
| <b>Molecular weight</b>        | 64.303 kDa  |
| <b>Expression System</b>       | HEK293  |
| <b>Species</b>                 | Human   |
| <b>Physical Appearance</b>     | Liquid  |
| <b>Buffer</b>                  | 1xPBS, pH7.4  |
| <b>Concentration</b>           | 0.800 mg/ml, BSA was used as standard and protein concentration was determined by Bradford                                      |
| <b>Stability &amp; Storage</b> | For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)  |



**Thank you for choosing Detai Biological. We are dedicated to serving you!**  
**For Research Use Only**