


Kazal-type serine protease inhibitor domain-containing protein 1 Protein

Cat.No: DTD0046

DESCRIPTION

Name	Kazal-type serine protease inhibitor domain-containing protein 1 protein
Purification	The soluble target protein Kazal-type serine protease inhibitor domain-containing protein 1 was derived from secreted expression and obtained by affinity chromatography purification
Purity	>85% by SDS-PAGE
Molecular weight	59 kDa
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1xPBS, 5% Glycerol, pH7.4
Concentration	1.400 mg/ml, BSA was used as the standard, and the protein concentration was determined by the Bradford method
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only