


Cerebellin-1 Protein

Cat.No:DTD0071

DESCRIPTION

Name	Cerebellin-1 protein
Purification	The soluble target protein Cerebellin-1 was derived from secreted expression and obtained by affinity chromatography purification
Protein Concentration	At 1.500 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA
Molecular weight	21 kDa
Purity	90% by SDS-PAGE
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
 For Research Use Only