


# Lysozyme C Protein

Cat.No:DTD0089

## DESCRIPTION

|                                |   |
|--------------------------------|---|
| <b>Name</b>                    | Lysozyme C protein  |
| <b>Purification</b>            | Lysozyme C protein was obtained from intracellular expression followed by purification steps. |
| <b>Solubility</b>              | Recommend to reconstitute the lyophilized LZ83 in sterilized water (< 160 µg/ml).             |
| <b>Molecular weight</b>        | 19 kDa  |
| <b>Purity</b>                  | >90% by SDS-PAGE  |
| <b>Species</b>                 | Chimpanzee  |
| <b>Expression System</b>       | HEK293  |
| <b>Buffer</b>                  | 1×PBS, pH7.4  |
| <b>Stability &amp; Storage</b> | For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)                |



**Thank you for choosing Detai Biological. We are dedicated to serving you!**  
**For Research Use Only**