


Butyrophilin-like protein 2

Cat.No:DTD0117

DESCRIPTION

Name	Butyrophilin-like protein 2
Purification	The soluble target protein Butyrophilin-like protein 2-His was derived from secretory expression and obtained by affinity chromatography purification
Purity	>90% by SDS-PAGE
Molecular weight	48 kDa
Expression System	HEK293
Species	Mouse
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Concentration	At 1.25 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA assay
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only