


Cobalamin binding intrinsic factor Protein

Cat.No:DTD0125

DESCRIPTION

Name	Cobalamin binding intrinsic factor protein
Purification	The soluble target protein Cobalamin binding intrinsic factor was derived from secreted expression and obtained by affinity chromatography purification
Protein Concentration	At 0.280 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA
Molecular weight	46 kDa
Purity	>90% by SDS-PAGE
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only