


Growth hormone receptor Protein

Cat.No:DTD0131

DESCRIPTION

Name	Growth hormone receptor protein
Purification	The soluble target protein Growth hormone receptor was derived from secreted expression and obtained by affinity chromatography purification
Protein Concentration	At 0.900 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA
Molecular weight	28 kDa
Purity	>85% by SDS-PAGE
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only