


Pancreatic secretory granule membrane major glycoprotein GP2

Cat.No:DTD0315

DESCRIPTION

Name	Pancreatic secretory granule membrane major glycoprotein GP2
Purification	The soluble target Pancreatic secretory granule membrane major glycoprotein GP2 was derived from secreted expression and was purified by affinity chromatography
Purity	>85% by SDS-PAGE
Molecular weight	54.4 kDa
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, 5% Glycerol, pH7.4
Concentration	0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only