


Green fluorescent protein

Cat.No:DTD0345

DESCRIPTION

Name	Green fluorescent protein
Purification	The soluble target protein Green fluorescent protein was derived from secreted expression and was purified by affinity chromatography
Purity	>85% by SDS-PAGE
Molecular weight	27.7 kDa
Expression System	HEK293
Species	Aequorea victoria
Physical Appearance	Liquid
Buffer	1×PBS, 5% Glycerol, pH7.4
Concentration	0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)

Restriction Digestion Map



Thank you for choosing Detai Biological. We are dedicated to serving you!
 For Research Use Only