


Nucleo protein

Cat.No:DTD0357

DESCRIPTION

Name	Nucleo protein
Purification	The soluble target Nucleo protein was derived from secreted expression and was purified by affinity chromatography
Purity	>85% by SDS-PAGE
Molecular weight	49 kDa
Expression System	HEK293
Species	Dobrava-Belgrade orthohantavirus
Physical Appearance	Liquid
Buffer	1×PBS, 5% Glycerol, pH7.4
Concentration	0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only