


Breast cancer type 1 susceptibility Protein

Cat.No:DTD0376

DESCRIPTION

Name	Breast cancer type 1 susceptibility protein
Purification	The soluble target protein Breast cancer type 1 susceptibility protein was derived from secreted expression and was purified by affinity chromatography
Purity	>85% by SDS-PAGE
Molecular weight	34.4 kDa
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, 5% Glycerol, pH7.4
Concentration	0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only