


RecQ-like DNA helicase BLM Protein

Cat.No:DTD0384

DESCRIPTION

| | |
|--------------------------------|--|
| Name | RecQ-like DNA helicase BLM protein |
| Purification | The soluble target protein RecQ-like DNA helicase BLM was derived from secreted expression and was purified by affinity chromatography |
| Purity | >85% by SDS-PAGE |
| Molecular weight | 75.2 kDa |
| Expression System | HEK293 |
| Species | Human |
| Physical Appearance | Liquid |
| Buffer | 1×PBS, 5% Glycerol, pH7.4 |
| Concentration | 0.500 mg/ml, BSA was used as standard and protein concentration was determined by Bradford |
| Stability & Storage | For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws) |

SDS-PAGE & Western Blot



Thank you for choosing Detai Biological. We are dedicated to serving you!

For Research Use Only