


Cyno BCMA Protein

Cat.No:DTD0631

DESCRIPTION


Name	Cyno BCMA protein
Purification	The soluble target protein Cyno-BCMA was derived from secreted expression and was purified by affinity chromatography
Purity	>95% by SDS-PAGE; >95% by HPLC
Molecular weight	30 kDa
Expression System	HEK293
Tag	Mouse Fc
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Concentration	4 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)

SDS-PAGE




Lane1: Cyno BCMA -mFc
Lane2: BSA
M:SDS-PAGE marker

HPLC



ELISA



Thank you for choosing Detai Biological. We are dedicated to serving you!

For Research Use Only