


Human IL2R β & IL2R γ Protein

Cat.No:DTD0633


DESCRIPTION

Name	Human IL2R β & IL2R γ Protein
Purification	The soluble target protein Human IL2R β & IL2R γ was derived from secreted expression and was purified by affinity chromatography
Purity	Human IL2R β & IL2R γ Heterodimer Protein, Fc Tag&Fc Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%; >90% by HPLC
Molecular weight	60 kDa; 80kDa
Expression System	HEK293
Tag	Mouse Fc
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Concentration	3 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)


SDS-PAGE



HPLC



ELISA



Thank you for choosing Detai Biological. We are dedicated to serving you!

For Research Use Only