


Human TNFSF4 Protein

Cat.No:DTD0639


DESCRIPTION

Name	Human TNFSF4 protein
Purification	The soluble target protein Human TNFSF4 was derived from secreted expression and was purified by affinity chromatography
Purity	>95% by SDS-PAGE
Molecular weight	20 kDa
Expression System	HEK293
Tag	His
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Concentration	1.2 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)

SDS-PAGE



FACS



Thank you for choosing Detai Biological. We are dedicated to serving you!

For Research Use Only