


Human PD-L1 Protein

Cat.No:DTD0640


DESCRIPTION

Name	Human PD-L1 protein
Purification	The soluble target protein Human PD-L1 was derived from secreted expression and was purified by affinity chromatography
Purity	>95% by SDS-PAGE; >95% by HPLC
Molecular weight	30 kDa
Expression System	HEK293
Tag	His
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Concentration	0.8 mg/ml, BSA was used as standard and protein concentration was determined by Bradford
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)


SDS-PAGE



HPLC



ELISA



Thank you for choosing Detai Biological. We are dedicated to serving you!

For Research Use Only