

Flag-Tag Rabbit Monoclonal Antibody

Cat.No:DTF02

DESCRIPTION

Product Name	Flag-Tag Rabbit Monoclonal Antibody
Host Species	Rabbit
Ig class	IgG
Conjugate	Unconjugated
Purity	>95% (purified by proteinA affinity column)
Applications	WB / ELISA / IF / IP
Specificity	DYKDDDDK Tag Antibody, mAb, Rabbit, recognizes C-terminal, N-terminal and internal tagged fusion proteins.
Immunogen	DYKDDDDK synthetic peptide conjugate to KLH
Storage Buffer	PBS, pH7.4
Storage	For up to 12 months at -20°C or lower, After receive, the antibody should be divided into the right amount. Avoid repeated freezing and thawing cycles
Recommended Dilutions	WB: 1 : 500-1 : 10000 ELISA: 1 : 1000-1 : 10000 IP: 1 : 200-1 : 2000

