

Rabbit Anti GST Tag Antibody, pAb

Cat.No: DTG0201

DESCRIPTION

| | |
|----------------------|--|
| Name | Rabbit Anti GST Tag Antibody, pAb |
| Immunogen | Recombinant Glutathione-S-Transferase(GST) |
| Host Species | Rabbit |
| Specificity | This antibody recognizes GST tag fusion proteins. |
| Form | Liquid |
| Concentration | 0.5 mg/ml |
| Storage Buffer | 1×PBS, pH7.4, containing 50% glycerol and 0.03% Proclin 300. |
| Purification | Antigen affinity purification |
| Clonality | Polyclonal |
| Isotype | Rabbit IgG |
| Shipping and Storage | Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles. |
| Expiration Date | Specified on product label. |

Applications

| Application | Recommended Usage |
|--------------------|-------------------|
| ELISA | 0.005-0.5μg/ml |
| Western Blot | 0.1-1μg/ml |
| Other applications | User-optimized |

Examples




Figure 1. Detection of GST tag in recombinant fusion protein were analyzed by Western blot using Rabbit Anti GST Tag Antibody, pAb(Cat.DTG0201, at 1 μ g/mL concentration).

The signal was developed with Goat Anti-Rabbit IgG(H&L)[HRP], pAb(Cat.DTA0252).

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!