

# Rabbit Anti GST Tag Antibody[HRP], pAb

Cat.No: DTG0206


## DESCRIPTION

Name	Rabbit Anti GST Tag Antibody[HRP], pAb
Immunogen	Recombinant Glutathione-S-Transferase(GST)
Conjugate	Horseradish peroxidase (HRP)
Host Species	Rabbit
Specificity	This antibody recognizes GST tag fusion proteins.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 50% glycerol and 0.03% Proclin 300.
Purification	Antigen affinity purification
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

## Applications

Application	Recommended Usage
ELISA	0.05-0.5µg/ml
Western Blot	0.2-1µg/ml
Other applications	User-optimized

## Examples



**Figure 1.** Detection of GST tag in recombinant fusion protein were analyzed by Western blot using Rabbit Anti GST Tag Antibody[HRP], pAb(Cat.DTG0206, at 1  $\mu$ g/mL concentration).

感谢您选择德泰生物，我们竭诚为您服务！

**For Research Use Only!**