

His-Tag Rabbit Monoclonal Antibody

Cat. No: DTH02

DESCRIPTION

ProductName	His-Tag Rabbit Monoclonal Antibody
Host Species	Rabbit
Ig class	IgG
Purity	>95% (purified by proteinA affinity column)
Applications	WB / ELISA / IF / IP
Specificity	His-Tag Rabbit Monoclonal antibody detects endogenous levels of internal, C-terminal, or N-terminal 6His-tagged proteins.
Storage Buffer	PBS, pH7.4
Storage	For up to 12 months at -20°C or lower, After receive, the antibody should be divided into the right amount. Avoid repeated freezing and thawing cycles
Recommended Dilutions	WB: 1 :2500-1 :5000 ELISA: 1 :2500-1 :10000 IF: 1 :500