

Rabbit Anti c-Myc Tag Antibody[HRP], pAb

Cat.No: DTM0205

DESCRIPTION

Name	Rabbit Anti c-Myc Tag Antibody[HRP], pAb
Immunogen	c-Myc epitope tag peptide conjugated to keyhole limpet haemocyanin(KLH)
Conjugate	Horseradish peroxidase (HRP)
Host Species	Rabbit
Specificity	This antibody recognizes C-terminal, N-terminal, and internal region of c-Myc tag fusion proteins.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 50% glycerol and 0.03% Proclin 300.
Purification	Antigen affinity purification
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	0.05-0.5 µg/ml
Western Blot	0.2-1µg/ml
Other applications	User-optimized

Examples




Figure 1. Detection of c-Myc tag in recombinant fusion protein were analyzed by Western blot using Rabbit Anti c-Myc Tag Antibody[HRP], pAb (Cat.DTM0205, at 1 μ g/mL concentration).

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!