

Rabbit Anti-Alpaca IgG [HRP], pAb

Cat:DTA0276

DESCRIPTION

Name	Rabbit Anti-Alpaca IgG [HRP], pAb
Immunogen	Alpaca IgG
Conjugate	Horseradish peroxidase (HRP)
Host Species	Rabbit
Specificity	This antibody recognizes Alpaca IgG, including IgG1, IgG2, IgG3. And also has good reactivity to llama and Camel IgG1, IgG2, IgG3.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 50% glycerol and 0.03% Proclin 300.
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:5,000-1:50,000
Western Blot	1:1,000-1:10,000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!