

Rabbit Anti-Alpaca IgG[iFluor 488],pAb

Cat:DTA0279

DESCRIPTION

Name	Rabbit Anti-Alpaca IgG[iFluor 488],pAb
Immunogen	Alpaca IgG
Conjugate	iFluor 488
Host Species	Rabbit
Specificity	This antibody recognizes Alpaca IgG, including IgG1, IgG2, IgG3. And also has good reactivity to Llama and Camel IgG1, IgG2, IgG3.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 50% glycerol and 0.03% Proclin 300.
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
Flow Cytometry	1:500-1:2,000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!