

Rabbit Anti-Camel IgG[FITC],pAb

Cat:DTA0289

DESCRIPTION

Name	Rabbit Anti-Camel IgG[FITC],pAb
Immunogen	Camel IgG
Conjugate	FITC
Host Species	Rabbit
Specificity	This antibody recognizes Camel IgG, including IgG1, IgG2, IgG3. And also has good reactivity to Alpaca and Llama IgG1, IgG2, IgG3.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 50% glycerol and 0.03% Proclin 300.
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
Flow Cytometry	1:500-1:2,000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!