

Rabbit Anti-Camelid VHH[HRP],pAb

Cat:DTA0291

DESCRIPTION

Name	Rabbit Anti-Camelid VHH[HRP],pAb
Immunogen	Camelid VHH
Conjugate	Horseradish peroxidase (HRP)
Host Species	Rabbit
Specificity	This antibody recognizes camelid VHH, recombinant VHH and humanized VHH.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 1% BSA 50% glycerol and 0.03% Proclin 300.
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:5,000–1:50,000
Western Blot	1:1,000–1:10,000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!