

# Rabbit Anti-Camelid VHH[iFluor 488],pAb

Cat:DTA0294

## DESCRIPTION

|                      |                                                                                                                                                              |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name                 | Rabbit Anti-Camelid VHH[iFluor 488],pAb                                                                                                                      |
| Immunogen            | Camelid VHH                                                                                                                                                  |
| Conjugate            | iFluor 488                                                                                                                                                   |
| Host Species         | Rabbit                                                                                                                                                       |
| Specificity          | This antibody recognizes camelid VHH, recombinant VHH and humanized VHH.                                                                                     |
| Form                 | Liquid                                                                                                                                                       |
| Concentration        | 0.5 mg/ml                                                                                                                                                    |
| Storage Buffer       | 1×PBS, pH7.4, containing 1% BSA 50% glycerol and 0.03% Proclin 300.                                                                                          |
| Clonality            | Polyclonal                                                                                                                                                   |
| Isotype              | Rabbit IgG                                                                                                                                                   |
| Shipping and Storage | Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles. |
| Expiration Date      | Specified on product label.                                                                                                                                  |

## Applications

| Application        | Recommended Usage |
|--------------------|-------------------|
| Flow Cytometry     | 1:500-1:2,000     |
| Other applications | User-optimized    |

感谢您选择德泰生物，我们竭诚为您服务！

**For Research Use Only!**